

# *DaVinci connects to iPad & iPhone*

DaVinci™ 

## The Revolutionary All-in-One HD CCTV

with Text-to-Speech(OCR) and 3 camera viewing positions

See near, far and everything in between.

- **Sony® Full HD 1080p camera** - with full HD you will experience amazing picture detail
- **Auto focus camera** - for reading, distance viewing or self-viewing
- **Reads aloud with Nuance® software** - read any printed text aloud with the push of a button
- **24" high resolution HD LCD** - for crystal clear pictures and vibrant colors
- **28 adjustable viewing modes allow you to personalize DaVinci** - perfect for reading magazines, newspapers, etc.
- **Magnification** - up to 77x



### Here is all you need:

Apple iPad Dock Connector to VGA (female) Part # MC552ZM/B or if you are using 4th generation iPad you need a cable called Lighting to VGA (female) MD825ZM/A. Also required is a VGA cable (male to male). One end connects to the Apple cable and the other connects to the VGA connector at the back of the DaVinci or Merlin elite monitor.



**CALL NOW FOR MORE INFORMATION AND YOUR FREE DEMONSTRATION**

# 888-811-3161

[www.enhancedvision.com](http://www.enhancedvision.com)